

PRO COMFORT

Inhaler Spacer Chamber



Details:

- An inhaler spacer chamber also known as an Anti-static valved holding chamber.
- Designed to provide optimal delivery of your inhaler medication.
- Comes in Child and Adult Sizes
- Compatible with most MDI (metered dose inhaler) actuators.
- Helps keep you from spraying drug on the back of your throat.

Size Details:

- Chamber Size:
5.7cm Width x 14.9cm Height
- Adult Mask size:
12.5cm mask Width x 8.8cm Height
- Child Mask size:
10cm Mask Width x 6.3cm Height
- Transparent mouthpiece helps coordinate timing of medication inhalation.
- Valve and end cap remove easily for cleaning, and the valve can be replaced to extend the life of your chamber.
- Helps eliminate unpleasant tastes of some medications.
- Low resistance valve for easy inhalation as well as exhalation.

| Item # | UPC # | NDC # | BRAND NAME | AWP/SWP | Box QTY. | Case QTY. |
|--------|--------------|---------------|---------------------------|---------|----------|-----------|
| 7287 | 850632007287 | 50632-0007-28 | Child Mask Spacer Chamber | \$83.71 | each | 50 |
| 7270 | 850632007270 | 50632-0007-27 | Adult Mask Spacer Chamber | \$83.71 | each | 50 |



Home Aide Diagnostics, Inc.

1072 S Powerline Rd., Deerfield Beach, FL 33442

www.homeaide.us | 1-800-915-0116