



# Knee Support

## Knee Brace

### BENEFITS & USES

- o Thermal neoprene sleeves which provide therapeutic body heat retention while providing firm compression and soft tissue support.
- o Hinges on both sides of the knee.
- o Adjustable Velcro® straps.
- o Fits right or left knee.
- o Rigid uprights for additional medial-lateral support.
- o Provides stabilization, support and compression for mild knee sprains and ligament injuries.
- o The ultimate brace for weak or injured knees.
- o For mild medial/lateral support of the knee.
- o Stabilize collateral ligaments, Protect meniscus, Decrease pain symptoms, Improve joint function, Stabilize knee joint.
- o Helps with Arthritic knee pain.

